



# eATC-S Automated Thermal Cycler

ULTRASSAY BIOTECH CO., LTD.

2022.07

# Index...

1. [Design concept...](#)
2. [Product appearance...](#)
3. [Exclusive technology...](#)
4. [Installation dimension...](#)
5. [System connection...](#)
6. [Communication protocol and host computer interface...](#)
7. [Repair kit...](#)
8. [Fast repair...](#)
9. [Interchangeability verification...](#)
10. [Main technical parameters...](#)
11. [Automated series products...](#)
12. [Pain points and solutions...](#)
13. [Competitive analysis...](#)

# 1. Design concept...






## eATC-S Automated thermal cycler




- **Design concept** : Focus on solving automated users' concerns about the reliability and maintainability of PCR instruments.
- **Product spec** : 96/384 optional blocks, apply to full skirt plate and short tube only.
- **Unique function** : Users can quickly change the PCR module independently, free of temperature calibration, leakage prevention.
- **Exclusive technology** : Rapid insertion and removal technology, free of calibration module technology; Leak-proof sealing bracket technology; liquid cooling technology.
- **Market positioning** : Used for supporting automatic workstations.



# 2. Product appearance...



# 3. Exclusive technology...




## Rapid insertion and removal technology

- This technology allows customers to replace the PCR module themselves under the guidance of the manufacturer.

**Free of calibration**

## Free of calibration module technology

- This technology solves the problem of temperature calibration after the replacement of thermal module, and makes it possible for users to maintain independently.



## Anti-evaporation sealing bracket technology

- The technology solves the leakage problem when the robot arm grabs the gasket after the amplification is completed, and effectively prevents the cross contamination of the amplification products.




**Liquid cooling technology**

## Liquid cooling remote heat dissipation technology


- It guides the heat of the thermal cycle process to the outside of the workstation for heat dissipation, effectively preventing the heat accumulation in the workstation and the pollution risk caused by air flow disturbance.
- It makes PCR hosts to be tightly arranged in workstations without worrying about air cooling blocking the vents.



# 4. Installation dimension...




Dimension of host (Thermal cycler unit)



# 5. System connection...

Connection diagram





# 6. Communication protocol and host interface...



Serial communication protocol manual

Automated PCR for workstation (eATC-S)

End running

Running log

Create program

Select serial port

End running

Ruuning.....



# 7. Repair kit ...



Pluggable thermal cycler module

Heat conducting silicone grease


hexagon wrench

Silicone coating roller

- Peltier is the core component of PCR instrument, but also damageable component, it determines the operation of PCR instrument.
- Replacing Peltier is a very professional job, which requires a lot of professional knowledge and special tools.
- This toolbox solves the above problems. We integrate the design of block, temperature sensor, circuit, plug, Peltier and pressing mechanism to form a standard thermal module, which is packed in a toolbox with thermal interface materials and disassembly tools for spare parts. When necessary, the user can complete the replacement in a few minutes without temperature calibration.
- The "Quick repair kit box" needs to be bought separately.



# 8. Fast repair...



# 9. Interchangeability verification...

1. Repeated disassembly and assembly on the same module and the same host...**Completely interchangeable**

2. Repeated disassembly and assembly on the same module and the different host.....**Completely interchangeable**

0	1	2	3	4	5
0	"参数"	"测试点"	"实测数据"	"单位"	"内控校准目标"
1	"最大升温速率"	"数据过程瞬间"	4.2	"°C/s"	"..."
2	"最大降温速率"	"数据过程瞬间"	7	"°C/s"	"..."
3	"平均升温速率"	"50-90°C"	1.68	"°C/s"	1.6
4	"平均降温速率"	"90-50°C"	1.65	"°C/s"	1.5
5	"温度准确度"	"95°C"	-0.02	"°C"	0.3
6	"温度准确度"	"72°C"	0.02	"°C"	0.3
7	"温度准确度"	"55°C"	0.03	"°C"	0.3
8	"模块温度均匀性"	"95°C"	0.46	"±/°C"	0.6
9	"模块温度均匀性"	"72°C"	0.32	"±/°C"	0.6
10	"模块温度均匀性"	"55°C"	0.17	"±/°C"	0.4
11	"模块温控精度"	"95°C"	0.02	"±/°C"	0.2
12	"模块温控精度"	"72°C"	0.01	"±/°C"	0.2
13	"模块温控精度"	"55°C"	0.05	"±/°C"	0.2
14	"低温冷藏(准确度)"	"4°C"	0.02	"°C"	1

0	1	2	3	4	5
0	"参数"	"测试点"	"实测数据"	"单位"	"内控校准目标"
1	Module TEMP uniformity	"95°C"	0.46	"±/°C"	0.6
2	Module TEMP uniformity	"77°C"	0.37	"±/°C"	0.6
3	Module TEMP uniformity	"55°C"	0.17	"±/°C"	0.4

0	1	2	3	4	5
0	"参数"	"测试点"	"实测数据"	"单位"	"内控校准目标"
9	"模块温度均匀性"	"72°C"	0.32	"±/°C"	0.6
10	"模块温度均匀性"	"55°C"	0.17	"±/°C"	0.4
11	"模块温控精度"	"95°C"	0.02	"±/°C"	0.2
12	"模块温控精度"	"72°C"	0.02	"±/°C"	0.2
13	"模块温控精度"	"55°C"	0.05	"±/°C"	0.2

0	1	2	3	4	5
0	Module TEMP uniformity	"95°C"	0.46	"±/°C"	0.6
1	Module TEMP uniformity	"77°C"	0.37	"±/°C"	0.6
2	Module TEMP uniformity	"55°C"	0.17	"±/°C"	0.4

0	1	2	3	4	5
0	"参数"	"测试点"	"实测数据"	"单位"	"内控校准目标"
1	"最大升温速率"	"数据过程瞬间"	3.1	"°C/s"	"..."
2	"最大降温速率"	"数据过程瞬间"	2.5	"°C/s"	"..."
3	"平均升温速率"	"50-90°C"	2.05	"°C/s"	1.6
4	"平均降温速率"	"90-50°C"	1.83	"°C/s"	1.5
5	"温度准确度"	"95°C"	0.79	"°C"	0.3
6	"温度准确度"	"72°C"	-0.56	"°C"	0.3
7	"温度准确度"	"55°C"	-0.65	"°C"	0.3
8	"模块温度均匀性"	"95°C"	0.42	"±/°C"	0.6
9	"模块温度均匀性"	"72°C"	0.25	"±/°C"	0.6
10	"模块温度均匀性"	"55°C"	0.19	"±/°C"	0.4
11	"模块温控精度"	"95°C"	0.02	"±/°C"	0.2
12	"模块温控精度"	"72°C"	0.03	"±/°C"	0.2
13	"模块温控精度"	"55°C"	0.03	"±/°C"	0.2
14	"低温冷藏(准确度)"	"4°C"	-1.09	"°C"	1

0	1	2	3	4	5
0	Module TEMP uniformity	"95°C"	0.46	"±/°C"	0.6
1	Module TEMP uniformity	"77°C"	0.37	"±/°C"	0.6
2	Module TEMP uniformity	"55°C"	0.17	"±/°C"	0.4

0	1	2	3	4	5
0	"参数"	"测试点"	"实测数据"	"单位"	"内控校准目标"
9	"模块温度均匀性"	"72°C"	0.28	"±/°C"	0.6
10	"模块温度均匀性"	"55°C"	0.18	"±/°C"	0.4
11	"模块温控精度"	"95°C"	0.03	"±/°C"	0.2
12	"模块温控精度"	"72°C"	0.03	"±/°C"	0.2
13	"模块温控精度"	"55°C"	0.03	"±/°C"	0.2

0	1	2	3	4	5
0	Module TEMP uniformity	"95°C"	0.46	"±/°C"	0.6
1	Module TEMP uniformity	"77°C"	0.37	"±/°C"	0.6
2	Module TEMP uniformity	"55°C"	0.17	"±/°C"	0.4

# 10. Main technical parameters...

## Functions: With all the functions of our ETC821 thermal cycler

- Adapted blocks: 0.2ml\*96 wells full skirt plate, 384 wells plate
- Temperature accuracy:  $\leq \pm 0.2^{\circ}\text{C}$
- Block temperature control accuracy:  $\leq \pm 0.1^{\circ}\text{C}$
- Block temperature uniformity:  $\leq \pm 0.3^{\circ}\text{C}$  @55°C
- Max. heating rate:  $\geq 4^{\circ}\text{C/s}$
- Avg. heating rate:  $\geq 2^{\circ}\text{C/s}$  (50-90°C)
- Avg. cooling rate:  $\geq 1.5^{\circ}\text{C/s}$  (90-50°C)
- Block temperature control technology: 3 zones Peltier temperature control
- Cooling technology: Liquid cycle cooling
- Hot lid: electric hot lid

## Unique technology:

- Liquid cooling technology, fast module technology, free calibration technology.

## Experiment management:

- Administrator, visitor, password management, screen lock, log export.

## Expenimental safety:

- Peltier intelligent life warning, module DIY, module free calibration, quick repair kit.
- Communication & Control:
- Communication interface: S232C protocol, DB9 mother type seat.
- Control system: With operating screen and software. Through communication protocol control, remote control of hot lid movement (opening and closing), coordination of running program (selecting program, running, stopping, etc.) and feedback of host status (host SN, program status, hot lid location, host fault code, motor fault, etc.) can be performed on PC.
- Reaction plate grasping mode: The reaction plate can be operated by the robot arm from the front, from the sides or from above.

## Power supply & Dimension:

- Power supply: 100-240V ~ , 50/60Hz, 1000VA
- Net weight: 14kg(host 6.9kg, cooling unit 7.1kg)
- Gross weight: about 17kg
- Carton size: L\*W\*H=540\*450\*385mm





# 11. Automated series products...



**eATC-S**  
Automated thermal  
cycler

- Only suitable for full skirt plate, small size, full functions.
- Suitable for placement in the workstation or multiple sets of placement.
- can be used as an independent PCR instrument.



**eATC**  
Automated thermal  
cycler

- Suitable for full size test tubes. Large volume, full function.
- Suitable for placement outside the workstation.
- can be used as an independent PCR instrument.



**ECU831**  
Cooling unit

- Suitable for general workstation, 4~70°C
- Communication interface.
- Modules can be customized.





# 12. Painpoints and solutions...



Major defects of automated PCR in the market



Ultrassay's solution



# 13. Competitive analysis...



## Ultrassay

### Features :

- Module replacement DIY , free of calibration.
- Water cooling, compact layout, no troubles of shielding the air outlet.
- Local control, strong anti-interference, will not cause communication interruption.
- Small size, about 3.5 blocks in size.
- With operating screen, can be used as an independent PCR instrument.
- Performace : uniformity $\leq$ 0.2°C.

## Inheco

### Features :

- Module replacement requires professionals.
- The air inlet in the case is easily shielded by the blocks, which affects the heat dissipation performance and temperature change speed of the instrument.
- Remote control, poor anti-interference, easy to interrupt the signal.
- The hot lid wire is easy to break.
- Small size, about 3 blocks in size.
- Without screen, cannot be used as an independent PCR instrument.
- Performace: uniformity $\leq$ 0.15°C.

## Thermo Fisher

### Features :

- Module replacement requires professionals.
- Air cooling, the air inlet in the case is easily shielded by the blocks, which affects the heat dissipation performance and temperature change speed of the instrument.
- Remote control, poor anti-interference, easy to interrupt the signal.
- Small size, about 3.5 blocks in size.
- Without screen, cannot be used as an independent PCR instrument.
- Performace: uniformity $\leq$ 0.2°C.







Accurate you can trust !



<http://ultrassay.com>    <https://www.ybotech.com>