



FA/FA-N SERIES Internal Calibration Analytical Balance

Function:

- ★ Internal auto calibration
- ★ LCD (White back light with black font)
- ★ One key/Time trigger internal calibration (FA-N ONLY)
- ★ Electromagnetic force sensor
- ★ Aluminum shell/Large space wind proof cover
- ★ Tare function/Counting/Percentage
- ★ Unit conversion(g/mg/ct/oz)...
- ★ Under weighing
- ★ Sensitivity/Speed set
- ★ Overload alarm/Fault alarm/Level indicator

FA/FA-N SERIES						
Model (External Calibration)	FA1004	FA1104	FA1204	FA2004	FA2104	FA2204
Model (Internal Calibration)	FA1004N	FA1104N	FA1204N	FA2004N	FA2104N	FA2204N
Capacity	100g	110g	120g	200g	210g	220g
Readability	0.0001g					
Min Weighing	0.0004g					
Repeatability	±0.0002g					
Linearity	±0.0003g					
Stable Time	≤3S					
Pan Size	Ø80mm					
Appearance Size	340x215x350mm					
Draft Shield Size	240x190x265mm					
Communication	RS232(Standard) / USB、RJ45 (Optional)					
Temperature	15-35℃					
Cal.Weight	100g			200g		
Option	Printer					
Packing Size	Inner box:415x285x430mm			Outer box:475x315x485mm		
Gross Weight	11kg					