




Option: static eliminator



Application:



FA/FA-N SERIES Internal Calibration Analytical Balance

Function:

- ★ Internal auto calibration
- ★ 5 Inch touch panel
- ★ One key/Time trigger internal calibration
- ★ Electromagnetic force sensor
- ★ Aluminum shell/Large space wind proof cover
- ★ Tare function/Counting/Density/Percentage
- ★ Unit conversion (g/mg/ct/oz)...
- ★ Under weighing
- ★ Sensitivity/Speed set
- ★ Overload alarm/Fault alarm/Level indicator
- ★ Option:Density kit (Solid/Liquid)

| FA-T SERIES | | | | | | |
|-------------------|--------------------------------|---------|---------|-------------------------|---------|---------|
| Model | FA1004T | FA1104T | FA1204T | FA2004T | FA2104T | FA2204T |
| Capacity | 100g | 110g | 120g | 200g | 210g | 220g |
| Readability | 0.0001g | | | | | |
| Min Weighing | 0.0004g | | | | | |
| Repeatability | ±0.0002g | | | | | |
| Linearity | ±0.0003g | | | | | |
| Stable Time | ≤3S | | | | | |
| Pan Size | Ø80mm | | | | | |
| Appearance Size | 340x215x350mm | | | | | |
| Draft Shield Size | 240x190x265mm | | | | | |
| Communication | RS232(Standard)/ USB(Optional) | | | | | |
| Temperature | 15-35℃ | | | | | |
| Calibration | Internal Calibration | | | | | |
| Cal.Weight | 100g | | | 200g | | |
| Option | Printer | | | | | |
| Packing Size | Inner box:415x285x430mm | | | Outer box:475x315x485mm | | |
| Gross Weight | 11kg | | | | | |