



Features

- Integrated with incubator and shaker to save space and cost.
- High quality steel shell, polished stainless steel chamber.
- Large LCD screen to display temperature and shaking speed.
- With operation data memory function to eliminate faulty operation.
- Non-volatile memory saves settings during a power outage and automatically restarts the unit as originally programmed after power is restored.
- Auto-stop operation when door is opened. Strong air spring rod with easy opening and closing.
- Brushless DC motor, more stable and reliable.
- With UV Sterilizer
- Equipped with leakage protection.

Specifications

Model	YB-HSI-211B	YB-HSI-211C
Shaking Speed (rpm)	20-300	
Speed Accuracy (rpm)	±1	
Swing Amplitude (mm)	Φ26	
Standard Configuration	500ml×28	2000ml×12
Maximum Capacity	250ml×36 or 500ml×28 or 1000ml×18	1000ml×18 or 2000ml×12 or 3000ml×8 or 5000ml×6
Tray Size (mm)	920×510	
Timing Range	1~9999min	
Temperature Range (°C)	4-60 (Cooling)	4-60(Cooling)
Temperature Accuracy (°C)	±0.1	
Temperature Uniformity (°C)	±1	
Display	LCD	
Tray Included	1	
Exterior Size(W×D×H)mm	120×74×80	120×74×100
Net Weight(kg)	174	183
Power Rating (W)	866	951
VolumeW×D×H(mm)	970×565×280 (155L)	970×565×480 (265L)
Power Supply	AC220V±10% , 50-60Hz	