



## JA-P SERIES Analytical Balance

Function:

- ★ LCD/AC and DC power supply
- ★ Tare function/Counting/Unit conversion(g/ct/oz)...
- ★ Min weighing set
- ★ Overload alarm/Level indicator
- ★ Option: Dual display/Interface/Printer/Dust cover

Application:



JA-P specifications				
Model	JA103P	JA203P	JA303P	JA103PL
Capacity	110g	210g	310g	110g/610g
Readability	0.001g			0.001g/0.01g
Min Weighing	0.001g-0.004g			1d-4d
Repeatability	±0.002g			±2d
Linearity	±0.003g			±3d
Stable Time	≤3S			
Pan Size	Ø90mm			Ø90mm/Ø110mm
Appearance Size	270x200x270mm			
Draft Shield Size	175x180x205mm			
Calibration	External Calibration			
Cal.Weight	100g	200g		
Packing Size	360x265x350mm/730x560x380mm			
Gross Weight	4kg/16kg			