

## VERTICAL CLASS B PRESSURE STEAM STERILIZER



### APPLICATION

MSL.B series vertical pressure steam sterilizer is suitable for sterilizing instruments, lumen instrument, textile, rubber, liquid, wastes and more, widely used in the clinic, hospital, laboratory and medical health center

### MAIN FEATURE

- Microcomputer control system, touch screen displaying
- Printer optional
- Built-in vacuum pump and vacuum drying system
- Built-in water tank, automatic sterilization process
- Built-in rapid steam generator
- Limit vacuum: -0.096Mpa
- Pressure displaying range: 0 ~ 0.4Mpa
- Temperature uniformity:  $\leq \pm 0.5^{\circ}\text{C}$
- Sterilization temperature range:  $105^{\circ}\text{C} \sim 138^{\circ}\text{C}$
- Chamber material: Stainless steel stainless steel
- Accessory: Stainless steel sterilizing basket

### MAIN TECHNICAL PARAMETER

Model	Chamber volume	Chamber dimensions (Φ×L)mm	Overall dimensions (L×W×H) mm	Net weight	Power supply
MSL.B50L	50L	316×667	480×620×990	80kg	3KW,AC 220V/110V
MSL.B80L	80L	386×695	550×690×1030	90kg	4.5KW,AC 380V/220V/110V
MSL.B100L	100L	386×860	550×690×1200	100kg	5.4KW,AC 380V/220V/110V

