

# MX-Mini

## Mini Vortex Mixer

### Compact fixed speed mixer



#### Features

- Compact, portable and fits nicely on bench-top or desk tops and table tops.
- Touch operation, fixed speed of 3000rpm for rapid mixing of samples.
- Powerful and efficient 12V DC brushless motor is smooth operation and has long life.
- Foot pad design secures the mixer in place and does not wobble during operation.



#### Specifications

Mixing motion	Orbital
Revolving diameter	4.8mm
Rotational speed	3000rpm( fixed)
Motor type	Brushless DC motor
Operating mode	Inching
Power supply	100~240V
Frequency	50/60HZ
Power	12W
Motor rating output	10W
Dimension[WxDxH]	100x100x79mm
Weight	0.5kg
Permissible ambient temperature and humidity	5~40°C, 80%RH
Protection class	IP21