

YB-2C SERIES

Precision Balance

- LCD/AC and DC power supply
- Tare function/Counting/Unit conversion (g/ct/oz)...
- Min weighing set
- Overload alarm/Level indicator
- Option: Dual display/Interface/Printer/ Windproof/Dust cover

Note: Capacity more than 1kg or accuracy less than 0.01g is not equipped with calibration weight



Application



Dual purpose



Gold Weighing



YB-1B/1BF SERIES

Precision Balance

- Fashional appearance
- Large size stainless steel pan (For YB-1BF)
- LCD (White back light with balk pot)
- AC and DC power supply
- Tare function/Counting/ Unit conversion (g/ct/oz)...
- Min weighing set
- Overload alarm/Level indicator
- Option: Dual display/Interface/Printer/ Dust cover



YB-1B



YB-1BF

YB-2C

Model	Capacity	Readability	Repeatability	Linearity	Pan Size	Stable Time
YB100-2C	110g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB200-2C	210g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB300-2C	310g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB500-2C	510g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB600-2C	610g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB1000-2C	1100g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB2000-2C	2100g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S
YB3000-2C	3100g	0.01g	±0.02g	±0.03g	Ø130mm	≤2S

YB-1B/1BF

Model	Capacity	Readability	Repeatability	Linearity	Pan Size	Stable Time
YB300-1B	310g	0.1g	±0.2g	±0.3g	Ø130mm	≤2S
YB600-1B	610g	0.1g	±0.2g	±0.3g	Ø130mm	≤2S
YB1000-1B	1100g	0.1g	±0.2g	±0.3g	Ø130mm	≤2S
YB2000-1B	2100g	0.1g	±0.2g	±0.3g	Ø130mm	≤2S
YB1000-1BF	1100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
YB2000-1BF	2100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
YB3000-1BF	3100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
YB5000-1BF	5100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
YB6000-1BF	6100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
YB10K-1BF	10000g	0.1g	±0.2g	±0.3g	160x160mm	≤2S