


YB-MX-T7 Series

Moisture Analyzer



Main Feature

- Electromagnetic Force Balance Sensor
- Accuracy 0.1mg
- Operator Authorization
- 200 datas storage
- One-key preheating
- High performance
- High standard
- Micro volume moisture analysis

Application



Medicine 5%



Chinese Medicine 7-15%



Plastic 0.25%



Stone Power 0.4-0.5%



Paper 6%



Sugar 5.5%

...

YB-MX-T7 Series				
Model	YB-1003MX-T7	YB-1003MX-T7*	YB-1004MX-T7 Pro	YB-1004MX-T7 Pro+
Capacity(g)	110	110	110	110
Accuracy (g)	0.001		0.0001	
Sensor	HBM		EFBS	
Display	7-inch touch screen			
Open Chamber	Manual	Automatic	Automatic	Automatic
Calibration	External	External	External	Internal
Min Weight (g)	0.004		0.0004	
Moisture Readability	0.01%		0.001%	
Temp. Range	40~200°C (230°C optional) 1°C step			
Stable Time	≤3s			
Pan Size	Φ96mm			
Operation Temp.	5-35°C			
Test Mode	Standard/Automatic/Mild/Customization			
Stop Mode	Automatic/Manual/Timming/Step/Aiming			
Heat Unit	Halogen Lamp			
Output	RS232、USB (Optional)			
Storage	200 datas, 200 test reports			
Other Functions	<ol style="list-style-type: none"> 1. Information input: Sample name/company/contact information 2. Adminstor/Operator Authorization/password 3. Time, date/ 200 groups data storage and review 4. Built-in test protocols 5. Various size of lables for printing 6. WIFI/Remote cloud services (Optional) 7. Chinese/English switch 8. GLP/GMP fromat storage 9. One-key preheating/calibration 10. Setting the cycle for automatic calibration (T7 Pro+ only) 11. Automatic door opening 12. Ultra-quite fan 			